

# INT STUDY RESTORATION

Shear bond strength of BRILLIANT Crios sandblasted with CoJet\*

R. Böhner, C. Kopfmann, A. Westendorf  
Internal Report Coltene 06/2016

## STUDY AIM


According instructions for use BRILLIANT Crios has to be sandblasted with  $\text{Al}_2\text{O}_3$  before application of ONE COAT 7 UNIVERSAL. The objective of this study was to show whether it is possible to use CoJet sandblasting system (3M ESPE) for pretreatment of BRILLIANT Crios before the application of ONE COAT 7 UNIVERSAL. A direct comparison of sandblasting with a commercial available  $\text{Al}_2\text{O}_3$  and CoJet should show whether the CoJet sandblasting system is appropriate for this purpose.

## EXPERIMENTAL SETUP

Samples of BRILLIANT Crios were either sandblasted with corundum (Renfert,  $\text{Al}_2\text{O}_3$ , 50  $\mu\text{m}$ ) or CoJet (3M ESPE, 30  $\mu\text{m}$ ). The sandblasted surface was bonded with ONE COAT 7 UNIVERSAL. SYNERGY D6 Flow was applied on the bonded surface. Bonding and Flow were cured through the BRILLIANT Crios material (3mm). Shear bond strength (SBS) was determined with a Watanabe setup after 24 h water storage at room temperature (n=8).


## RESULT

Both sandblasting materials, Renfert and 3M ESPE, show an edged surface. In both cases the BRILLIANT Crios surface is roughened after sandblasting. Shear bond strength of CoJet sandblasted samples are in the same range as the samples treated with Renfert  $\text{Al}_2\text{O}_3$ . In both cases the failure of the interface BRILLIANT Crios, ONE COAT 7 UNIVERSAL, SYNERGY D6 Flow was cohesive. Due the findings the CoJet sandblasting system (30  $\mu\text{m}$ , 3M ESPE) is an appropriate sandblasting material for pretreatment of BRILLIANT Crios before bonding with ONE COAT 7 UNIVERSAL.




\*CoJet is not a registered trademark of COLTENE


## BRILLIANT Crios



SEM pictures of Renfert and CoJet sandblasting material



SEM pictures of BRILLIANT Crios surfaces sandblasted with Renfert and CoJet



Shear bond strength (SBS) of bonded BRILLIANT Crios (ONE COAT 7 UNIVERSAL, SYNERGY D6 Flow) sandblasted with Renfert and CoJet