



Cotton Products
by ROEKO

Lasting absorbency


- highest quality
- smart solutions
- user friendly



Cotton Products

by ROEKO – a name – years of experience – a wealth of technical

Cotton products are used day in day out in the dental office and little thought is given to them. However the ROEKO products are renowned for their choice materials and excellent finish. Quality which pays off in the dental



Luna Gold

you see and feel the difference

Luna Gold has a mucophile cover which prevents contact between the fibres and the soft tissue. *Luna Gold* can therefore be removed without wetting the dental roll beforehand. There is no risk of tissue inflammation, an important feature when treating sensitive patients.

The exceptional quality of *Luna Gold* is easily recognised – the gold stripe differentiates *Luna Gold* from standard dental rolls.

For regular treatment we recommend

Luna – long lasting absorbency and stability

Cotton-roll – is extra soft and absorbs moisture quickly

Dental roll dispensers

ROEKO *Dental roll dispensers* keep your rolls within reach, while protecting them from spray and dust.

Roll-o-mat

The *Roll-o-mat* dispenser is simply ingenious: slight pressure on the lid is sufficient to dispense just one roll which stays on the fork ready for use.



Luna for kids

specially for very young patients

Standard dental rolls are often difficult to place in a child's mouth and can cause pain. *The Luna for kids* rolls are smaller and pleasanter for children, but like the standard *Luna* they keep their shape throughout treatment. They create sufficient working space even in the smallest mouth and their high absorbency prevents regular changing. Preventive dentistry, and treatment with fillings and sealants run more smoothly.

ROEKO dental rolls are non-chlorine bleached and packed in plastic bags protecting them from moisture and dust.

Lunamat

The *Lunamat* is equipped with a drawer. On opening the drawer 2-3 dental rolls are ready for use.



Parotisroll

multipurpose with spring action

Parotisroll is a unique dental roll with a spring action inside. Due to the spring action *Parotisroll* stays in place and provides more working space. *Parotisroll* exerts gentle pressure on the parotis gland reducing saliva flow. Used buccally, lingually or labially the spring action helps to keep the mouth open, the roll retracts the cheek, lips and tongue providing sufficient space for efficient and undisturbed treatment.


Parotisroll anatomic

by Dr. Hannelore Oberschachtsiek

Parotisroll anatomic is the *Parotisroll* with flattened ends. The flat ends make the placement of the dental roll easier and reduce the "foreign body feeling" for the patient. *Parotisroll anatomic* exerts less pressure onto the tongue's margin than standard dental rolls. The treatment of sensitive patients with a strong gag effect will be much easier.


Parotimat

Parotisroll dispenser holds 50 *Parotisroll*. The dispenser enables individual removal and hygienic storage.



al know how – 100% quality

production of high quality cotton products and their accessories demands a high level of technical know how. ental office, ensuring efficient and smooth procedures. Dive into ROEKO's incredible world of cotton.




Etching Pellets

ROEKO *Etching Pellets* are firmly pressed, cylindrically shaped pellets ideal for use with etching gels and pastes.

Cotton pellets

ROEKO *Cotton pellets* are super soft and extremely absorbent. They are ideal for drying cavities and applying medication. Fully automated production and best quality cotton ensure uniform size and shape from pellet to pellet.



Comprecap

Comprecap compression caps are a smart aid for the impression preparation. They stop bleeding naturally – by compression and control moisture.

Comprecap anatomic

Comprecap anatomic are anatomically formed compression caps with semicircular shaped spaces on two opposite sides and can therefore be easily placed on adjacent teeth.

Cotton pellet dispensers

ROEKO *Cotton pellet dispensers* keep your pellets within reach, while protecting them from spray and dust. Pellets can be removed with one hand, thanks to the stable design and easy-to-open lid. The rounded bars in ROEKO's pellet dispensers ensure the pellets keep their shape on removal, whereas wire mesh surfaces can destroy the pellets.

Solomat

an impressive eye catcher in your office and incredibly practical too

The spring base forces the cotton pellets to the surface. The rounded bars as opposed to wire mesh make it possible for the pellets to keep their shape on removal.

Solomat-N

for easy removal of single cotton pellets, etching pellets or foam plastic pellets

Autoclavable plastic dispenser with practical slot-in refill system. The cotton pellets are pushed to the surface by hand.

Solomat-Click

the smart pellet dispenser

Solomat-Click consists of grille, lid and the refill body. These three parts are simply clicked together when repurchasing the refill. The stainless steel base ensures stability.






Order Information

	Content	REF	
Luna	300 g		
	# 1	100 001	
	# 2	100 002	
	# 3	100 003	
	# 4	100 004	
	500 pcs.		
	# 1	105 001	
	# 2	105 002	
	# 3	105 003	
	1000 pcs.		
	# 1	104 001	
	# 2	104 002	
	# 3	104 003	
Luna for kids	200 pcs.		
# 2	100 032		
Luna Gold	300g		
	# 1	109 121	
	# 2	109 122	
	# 3	109 123	
	1000 pcs.		
	# 1	109 101	
	# 2	109 102	
# 3	109 103		
Cotton-roll	1000 pcs.		
# 1	104 011		
# 2	104 012		
# 3	104 013		
Roll-o-mat Plastic, white	1 pc.	120 004	
	134 x 100 x 48 mm		
Chrome-Nickel-Steel		120 005	
134 x 100 x 43,5 mm			
Can be used for dental roll sizes 1-3, standard length (36 - 40 mm). Not suitable for <i>Luna for kids</i> .			
Lunamat Plastic, white	1 pc.	120 000	
	52 x 106 x 100 mm		
Can be used for dental roll sizes 1-4, standard length (36 - 40 mm) and for <i>Luna for kids</i> .			
Parotisroll	∅	L	
	100 pcs.		
	# 1	9 mm	80 mm
	# 2	10 mm	80 mm
	# 3	10 mm	100 mm
	# 5	8 mm	100 mm
	anatomic	10 mm	100 mm
		50 pcs.	
# 160 001			
# 160 002			
# 160 003			
# 160 005			
# 160 013			
Parotimat	1 pc.	162 100	
11,2 x 10,3 x 8,2 cm			
Cotton pellets	∅		
	# 000	3,0 mm	1,5 g
	# 00	4,0 mm	4,0 g
	# 0	4,8 mm	4,0 g
	# 1	6,3 mm	10,0 g
	# 2	7,6 mm	10,0 g
			210 000
		210 100	
		210 110	
		210 111	
		210 112	

www.coltenewhaledent.com
© 2010 Coltène/Whaledent



Coltène/Whaledent Ltd.
The President Suite-A
Kendal House, Burgess Hill
Victoria Way
West Sussex, RH15 9NF/U.K.
Tel. +44 (0)1444 235486
Fax +44 (0)1444 870640
info@coltenewhaledent.co.uk

∅	Content	REF
Etching Pellets		
Small	2 mm	10,0 g
Medium	3 mm	10,0 g
Large	4 mm	10,0 g
		230 001
		230 002
		230 003
Solomat		
# 000	1 pc.	221 000
# 00		221 100
# 0		221 110
# 1		221 111
# 2		221 112
Solomat comes <u>without</u> cotton pellets		
Solomat-Click		
# 000	1 pc.	224 000
# 00		224 100
# 0		224 110
# 1		224 111
Solomat-Click comes <u>with</u> cotton pellets		
Refills		
# 000	3 pcs.	214 000
# 00		214 100
# 00		214 110
# 1		214 111
Solomat-N		
# 000	1 pc.	223 000
# 00		223 100
# 0		223 110
# 1		223 111
# 2		223 112
Solomat-N comes <u>with</u> cotton pellets		
Refills		
# 000	6 pcs.	213 000
# 00		213 100
# 0		213 110
# 1		213 111
# 2		213 112
Comprecap		
#1	120 pcs.	530 008
#1	600 pcs.	531 008
#3	120 pcs.	530 010
#3	600 pcs.	531 010
#5	60 pcs.	530 014
#5	400 pcs.	531 014
Set (#1/3/5)	150 pcs.	530 000
Comprecap anatomic		
#1	120 pcs.	531 201
#3	120 pcs.	531 203
#5	60 pcs.	531 205
Set (#1/3/5)	120 pcs.	531 299

Coltène/Whaledent AG
Feldwiesenstrasse 20
9450 Altstätten/Switzerland
Tel. +41 (0)71 757 53 00
Fax +41 (0)71 757 53 01
info@coltenewhaledent.ch

Coltène/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, Ohio 44223/USA
Tel. +1 330 916 8800
Fax +1 330 916 7077
info@coltenewhaledent.com

Coltène/Whaledent GmbH + Co. KG
Raiffeisenstrasse 30
89129 Langenau/Germany
Tel. +49 (0)7345 805 0
Fax +49 (0)7345 805 201
info@coltenewhaledent.de

coltène
whaledent